



## Abdominal Binders - Instructions

*Please read carefully.*

**THIS PRODUCT IS INTENDED FOR APPLICATION BY LICENSED HEALTH CARE PRACTITIONERS AS DIRECTED BY A PHYSICIAN OR OTHER MEDICAL AUTHORITY. PROPER FIT AND APPLICATION OF THIS PRODUCT MUST BE PERFORMED BY A MEDICAL PROFESSIONAL. THESE DIRECTIONS ARE GUIDELINES ONLY AND ARE NOT OFFERED AS MEDICAL RECOMMENDATIONS.**

### **Applying Abdominal Binders:**

Wrap binder around torso. Stretch until compression is introduced. Stick Velcro. For maximum compression apply in supine position.

### **Washing Instructions:**

Hand wash support in mild soapy cool water. Rinse thoroughly and Drip Dry on flat surface. **Do Not** use solvents or chlorine bleach. **Do Not** put in the dryer. **Do Not** iron.

### **Note:**

Always consult with a Health Care Practitioner before making any changes in your support.

*Discontinue use at once if any sign of dermatitis or other problem occurs. If discomfort occurs while wearing support, discontinue use at once and contact fitter or physician. Discontinue use at once if signs of prolonged wear affect usefulness of support. Contact fitter or physician with any questions about proper care or use. This item is intended for use in medical treatment, and not to be worn as a protective device during sports activity. **Support Should Not** be worn while sleeping. **Support Should Not** be worn by persons with phlebitis or circulatory problems. **Never** wear support over an open wound.*

