

Knee Orthoses with multi-position hinges



SIDE HINGE

Single axis hinge with adjustable extension stop at 0° - 15° - 30° - 45° to full ROM

| SIZE | KNEE CIRCUMFERENCE |
|---------------|--------------------|
| Ex-Small..... | 12" |
| Small..... | 13" |
| Medium..... | 14" |
| Large..... | 15" |
| Ex-Large..... | 16" |
| 2X-Large..... | 17" |
| 3X-Large..... | 18" |

MODEL K7 PDAC APPROVED L1832

ROM Hinged Knee Support

13" pull-on knee support with multi-position hinges provide control and support for ligament injuries and instabilities. 2 Velcro circumference straps with elastic compression and adjustable ROM extension control side hinges.

Adjustable ROM design features adjustable hinges for dial-in adjustment of extension control with stops at 0° - 15° - 30° - 45° to full ROM.