

# Patella Stabilizer With Hinge Options



## MODEL KJ-PS1 L1832

- 13" neoprene knee sleeve with open patella
- SIDE HINGE: Single axis hinge with adjustable extension stop at 0°- 15°- 30°- 45° to full ROM
- Exterior pressure strap and J-bar buttress
- Specify Left or Right



## MODEL KJ-PS2 L1820

- 90° back stop hinges
- 13" neoprene knee sleeve with open patella
- Exterior pressure strap and J-bar buttress
- Specify Left or Right



SIZE	KNEE CIRCUMFERENCE
Ex-Small.....	12"
Small.....	13"
Medium.....	14"
Large .....	15"
Ex-Large .....	16"
2X-Large.....	17"
3X-Large.....	18"