

# Alto™ Knee Orthoses

U.S. Orthotics Alto™ Knee Orthoses provide excellent control and support for ligament injuries and instabilities. An ideal brace for post MCL and ACL instabilities as well as medial and lateral instabilities. Our adjustable ROM design features polycentric adjustable hinges for convenient dial-in adjustment of flexion/extension control. Stops at 0°-15°- 30°- 60°- 90° to full ROM. Horizontal latch-top cap cover design ensures there is no inadvertent displacement of hinge pins. PDAC verified L1832



## PULL ON MODEL K901

13" length with easy pull-up and pull-off loops.

Flexion/extension hinge with contoured condyle pads.

2 Velcro circumference straps with elastic compression section.

### SIZE CHART

KNEE CENTER	
Ex-Small.....12"	Ex-Large.....16"
Small .....13"	2X-Large .....17"
Medium .....14"	3X-Large .....18"
Large.....15"	4X-Large.....19"



**ADJUSTABLE HINGE JOINT**  
Adjustable polycentric hinge joint for flexion/extension control with stops at 0°-15°- 30°- 60°- 90° positions to full ROM.

## WRAP-AROUND MODEL K-905

Anterior open wrap-around style for patients who cannot pull on a traditional pull-up support.

13" length with easy pull-up and pull-off loops.

Flexion/extension hinge with contoured condyle pads.

2 Velcro circumference straps with elastic compression section. Open popliteal.

### SIZE CHART

KNEE CENTER	
Small .....11"-13"	2X-Large .....19"-21"
Medium .....13"-15"	3X-Large .....21"-23"
Large.....15"-17"	4X-Large.....23"-25"
Ex-Large.....17"-19"	

 **U.S. Orthotics**® INC.

800-825-5228

www.usorthotics.com  
info@usorthotics.com

8605 Palm River Rd.

Tampa, FL 33619

Ph: 813-621-7797

Fax: 813-623-1055

Toll-Free: 800-825-5228