

Below Knee Neoprene Suspension Sleeves



ZIPPERED SLEEVE

ZIP SLEEVE™

Our SZ contoured sleeve comes with a side zipper for ease of application – excellent for the geriatric patient who can not pull on a traditional sleeve. 4" bottom is textured rubber inside for gripping prosthesis and top is nylon inside for comfort against the skin. 8" zipper with pull tab. Available with or without top strap. 100% neoprene rubber. Available in Black, Dark Brown, Caucasian. 4mm. 14" length. Specify left or right. Custom sizes available.



FLEXION SLEEVE SERIES

Our flexion sleeves have a natural 45 degree anatomical curve built in. This anatomical slope is more natural and will reduce bunching and extend the life of the sleeves. Either nylon or textured rubber inside. 5mm. 16 1/2" length. 100% neoprene rubber. Available in Black, Dark Brown, Caucasian. Custom sizes available.
SS45 Anatomically sloped sleeves with outside nylon and all textured rubber inside.
SS45Y Anatomically sloped sleeves with nylon inside and out.



STRAIGHT TUBULAR SLEEVES

Our tubular sleeves are basic tapered sleeves with nylon outside and either rubber or nylon inside. 5mm. 14" length. 100% neoprene rubber. Available in Black, Dark Brown, Caucasian. Custom sizes available.
SST Tubular sleeves all textured rubber inside.
SSTY Tubular sleeves all nylon inside.



S1 Nylon on the outside, textured rubber on the inside. Contour patella. 4mm. 14" length.

S2 Nylon inside and out. Contour patella. 4mm. 14" length. Custom sizes available.



S3 Nylon reinforced center section with textured rubber top and bottom sections. Contour patella. 4mm. 14" length. Custom sizes available.



S4 Nylon top section and a textured rubber bottom section with Velcro and elastic top strap for additional pressure. Contour patella. 4mm. 14" length. Custom sizes available.

 **U.S. Orthotics**[®] INC.

800-825-5228

www.usorthotics.com

orders@usorthotics.com

info@usorthotics.com

8605 Palm River Rd.
Tampa, FL 33619
Ph: 813-621-7797
Fax: 813-623-1055
Toll-Free: 800-825-5228

Night Splints
Elastic Back Supports
Neoprene Supports
Knee Immobilizers
Walking Boots
FULLBACK™ XPA

AFO
Maternity Belts
Ankle Supports
Wrist Supports
Side-Lace Corsets
Embroidered Lab Coats

Clavicle Straps
Arm Slings
Cervical Collars
Rib Belts & Binders
Prosthetic Suspension
ZIP-SLEEVE™

Pediatric Supports
SPINA™ IV LSO
SPINA™ II Plastic Orthoses
PACO™ Hyperextension Orthoses
Pre-Form Foot Orthotics
Shoulder Immobilizers