

# DiamondBack<sup>®</sup> LS Supports



## MODEL 85-Di L0627

Our Diamondback<sup>™</sup> lumbosacral support is back with anterior and posterior pre-molded inserts. Inserts can be re-molded. Easy hand-pull tabs on all ends. Additional top tension strap.

- 85-Di-8 8" back tapers to a 4" front
- 85-Di-10 10" back tapers to a 6" front
- 85-Di-12 12" back tapers to a 6" front

Sizes	Waist	Sizes	Waist
X-Small.....	25" – 29"	Ex-Large.....	41" – 45"
Small.....	29" – 33"	2X-Large.....	45" – 49"
Medium.....	33" – 37"	3X-Large.....	49" – 53"
Large.....	37" – 41"	4X-Large.....	53" – 57"

Customs upon request. Made in the U.S.A.