

# Knee Immobilizers



**MODEL MBGL L1830**  
Soft foam knee immobilizer.  
All wide alternating elastic straps.  
Side stays are adjustable with velcro for proper fit. Trimable foam.  
Universal circumference.  
Fits up to 27" thigh.

**7 Heights Available**

- MBGL-14
- MBGL-16
- MBGL-18
- MBGL-20
- MBGL-22
- MBGL-24
- MBGL-26



**MODEL MBE L1830**  
Deluxe universal circumference knee immobilizer with elastic patella cover.  
Foam can be trimmed.  
Alternating Velcro closures.  
Fits up to 29" thigh.

**4 Heights Available**

- MBE-16
- MBE-19
- MBE-22
- MBE-25

**MODEL MBE-W**

Wide MBE model available for up to 35" thigh.



**MODEL MBTR L1830**  
Tri-panel universal circumference knee immobilizer. Three panels allow for greater adjustability for patients with wide calf to thigh ratios. Fits up to 29" thigh.

**3 Heights Available**

- MBTR-16
- MBTR-19
- MBTR-22

MBTR Tri-Panel

