

# PACO™

## Plastic Anterior Control Orthosis

### SADMERC Verified L0472

A 3-point pressure orthosis that encourages a mild hyper-extended posture. Durable plastic lined with plush soft-foam. Double comfort padding at the sternal, back and lower abdominal areas – ideal for geriatric patients.

Height is measured from sternal notch to pubic – minus 5". Orthosis size is determined by brace height. Circumference is adjustable by straps at the back pad and ranges from 30" to 44".

#### Size Front Height

Ex-Small.....14" height  
 Small.....16" height  
 Medium.....18" height  
 Large.....20" height  
 Ex-Large.....22" height



LIGHTWEIGHT PLASTIC LINED WITH PLUSH SOFT-FOAM



3-POINT PRESSURE ORTHOSIS



DOUBLE COMFORT PADDING

 **U.S. Orthotics**® INC.

800-825-5228

[www.usorthotics.com](http://www.usorthotics.com)

8605 Palm River Rd.

Tampa, FL 33619

Ph: 813-621-7797

Fax: 813-623-1055

Toll-Free: 800-825-5228

[orders@usorthotics.com](mailto:orders@usorthotics.com)

Night Splints

Elastic Back Supports

Neoprene Supports

Knee Immobilizers

Walking Boots

FULLBACK™ XPA

[info@usorthotics.com](mailto:info@usorthotics.com)

AFO

Maternity Belts

Ankle Supports

Wrist Supports

Side-Lace Corsets

Embroidered Lab Coats

Clavicle Straps

Arm Slings

Cervical Collars

Rib Belts & Binders

Prosthetic Suspension

ZIP-SLEEVE™

Pediatric Supports

SPINA™ IV LSO

SPINA™ II Plastic Orthoses

PACO™ Hyperextension Orthoses

Pre-Form Foot Orthotics

Shoulder Immobilizers