

## OSB Ostomy belts

Elastic with reinforced stoma opening. Provides additional support to the abdomen around the stoma.

Choose color elastic ..... Black or white

Comes in 5 heights ..... 4" - 6" - 8" - 10" - 12

2 Sizes ..... Short 28" – 38"      Long 38" – 48" (Custom Upon request)


4 Access Hole Sizes 2 ½" – 3" – 3 ½" – 4"

BELTS come standard with a 3 centered opening.

All belts come with a 1/16" thin polypropylene plastic reinforced stoma opening.

### Standard stoma access hole placement

Standard design below shows 6" belt with 3" access hole opening placed in the center of garment.



### Custom stoma access hole placement

Select 4 steps.

- 1) Select belt height & color.
- 2) Select patient left or right.
- 3) Select access hole size.
- 4) Access hole location. From bottom edge to access hole bottom edge.

