

ALTO[®] ROM Knee Orthoses

PDAC Approved
 L1832 / L1833



MODEL K-964

Alto™ ROM Knee Orthoses – Wrap-Around with covered hinges

ALTO ROM Knee Orthoses Model K-906 is a 13" anterior open wrap-around style for patients who can not pull on a traditional pull up support.

2 circumference straps and flexion / extension side hinges provide excellent control and support for ligament injuries and instabilities.

Covered latch-top cap cover design ensures there is no inadvertent displacement of hinge pins.

An ideal brace for post MCL and ACL instabilities.

Our adjustable R.O.M. design features polycentric adjustable hinges for convenient dial-in adjustment of flexion/extension control.

ROM Stops at 0° – 15° – 30° – 60° – 90° to Full ROM



Sizes	Knee Circumference
Small.....	11" – 13"
Medium	13" – 15"
Large	15" – 17"
Ex-Large.....	17" – 19"
2X-Large.....	19" – 21"
3X-Large.....	21" – 23"
4X-Large.....	23" – 25"
Customs upon request.	