

ALTO™ Uni - Adjustable Knee Orthosis



K-908 Uni - 13" Height - 2 Straps



K-909 Uni - 17" Height - 4 Straps

MODEL K-908-U PDAC L1832
MODEL K-909-U PDAC L1832

Universal 3-Panel Adjustable ROM Hinged Knee Support

- * 13" and 17" wrap-around knee supports with multi-position hinges provide control and support for ligament injuries and instabilities.
- * Adjustable Velcro circumference straps.
- * Adjustable hinge pocket assemblies for proper medial/lateral placement.
- * ROM design features adjustable hinges for dial-in adjustment of extension control. Stops at 0° – 15° – 30° – 45° to full ROM.



Adjustable, trimmable panels to fit a wide range of thigh-to-calf ratios.



Sizing

Adjustable to fit range of:

Thigh	13" - 30"
Knee	11.5" - 28"
Calf	12" - 28"