

# V-Line<sup>®</sup> Spinal Orthoses Pro-V28 L0464



Our professional grade V-Line<sup>®</sup> Spinal Orthoses series features comfort, stability and a wide array of options. A single piece design makes for easy application and effective compression. Pre-molded posterior inserts come in 0° – 7° – 12° – 25° lordosis with a spinal relief cut-out. Pre-molded anterior inserts comes neutral or pendulous and are pre-slotted to accommodate an anterior thoracic extension. The anterior pocket is fully adjustable providing correct positioning and accommodating a wider range of sizes. Designed with two smooth cross-pull compression straps, with elastic tension ends. Hand-pull tabs on all ends.



U.S. Orthotics PRO-V Model V28 (L0464) spinal orthosis is a single piece TLSO design that accommodates appropriate motion restriction throughout the healing process and maximum control for your patient with severe, high-level back pain.

The Pro-V28 delivers triplanar control by means of a posterior-lateral panel and circumferentially-overlapping plastic, with an anterior panel; sternal plate assembly and padded shoulder straps.

Two smooth cross-pull compression straps with elastic tension ends make for easy application and effective compression.

Optional under the axilla padded straps or patient friendly over-the-shoulder straps that apply easy as a backpack.  
 Hand-pull tabs on all ends.

**Pro-V28**

**L0464**



SIZE	HIPS
Small	32" – 38"
Medium	36" – 42"
Large	40" – 46"
Ex-Lg	44" – 50"
2XL	48" – 54"
3XL	52" – 58"
4XL	56" – 62"