

Arm Slings



MODEL S110

ENVELOPE SLING WITH BUCKLE CLOSURE

Lightweight cloth with a long sleeve encasing the entire forearm and hand to help prevent ulna deviation of the hand. Three bar suspender buckle adjustment to fit any size patient. Metal hook and catch closure. For maximum comfort the webbing is padded at the neck. Fits both left and right.



MODEL S110V

ENVELOPE SLING WITH VELCRO CLOSURE

Lightweight cloth with a long sleeve encasing the entire forearm and hand to help prevent ulna deviation of the hand. Velcro® closure. For maximum comfort the webbing is padded at the neck. Fits both left and right.

Size/Both Models Envelope Dimensions

Universal	7 1/2" x 16 1/2"
Small.....	7" X 12"
Medium.....	8" X15"
Large.....	9" X 18"