

ALTO[®] ROM Knee Orthoses



MODEL K-906 L1832

Alto™ ROM Knee Orthoses – Wrap-Around

Alto™ ROM Knee Orthoses Model K-906 is a 13" anterior open wrap-around style for patients who can not pull on a traditional pull up support.

2 circumference straps and flexion/extension control side hinges provide excellent control and support for ligament injuries and instabilities.

Latch-top cap cover design ensures there is no inadvertent displacement of hinge pins. An ideal brace for post MCL and ACL instabilities.

Our adjustable R.O.M. design features polycentric adjustable hinges for convenient dial-in adjustment of flexion/extension control.

ROM Stops at 0° – 15° – 30° – 60° – 90° to Full ROM

NEOPRENE or DRY FIT Fabric

Sizes

| | |
|----------------|-----------|
| Small..... | 11" – 13" |
| Medium | 13" – 15" |
| Large | 15" – 17" |
| Ex-Large | 17" – 19" |
| 2X-Large..... | 19" – 21" |
| 3X-Large..... | 21" – 23" |
| 4X-Large..... | 23" – 25" |

Customs upon request.

Knee Circumference

Indications

Post-operative management & rehabilitation. Locked or limited motion control of knee. Mild to moderate ACL, PCL, MCL, LCL. Combined ligament instabilities. Medial and lateral collateral ligament instabilities.

