

# Anterior Opening ROM Hinged Knee Support



**SIDE HINGE**

Single axis hinge with adjustable extension stop at 0° - 15° - 30° - 45° to full ROM

| SIZE          | KNEE CIRCUMFERENCE |
|---------------|--------------------|
| Small.....    | 11"-13"            |
| Medium.....   | 13"-15"            |
| Large.....    | 15"-17"            |
| Ex-Large..... | 17"-19"            |
| 2X-Large..... | 19"-21"            |
| 3X-Large..... | 21"-23"            |
| 4X-Large..... | 23"-25"            |

**MODEL KW8 PDAC APPROVED L1832**

**Anterior opening wrap-around knee orthosis with multi-position hinges**

12" wrap-around knee support with multi-position hinges provide control and support for ligament injuries and instabilities. 2 Velcro circumference straps with elastic compression and adjustable extension control side hinges.

Adjustable ROM design features adjustable hinges for dial-in adjustment of extension control with stops at 0° – 15° – 30° – 45° to full ROM.