

# Lumbosacral Support Model 80



## LUMBOSACRAL MODEL 80


### L0625

Deluxe back support with a cushioned fabric back pad, elastic side panels and Velcro closures. Additional 4" elastic top tension strap. 2 enclosed side stays help prevent roll-over and 4 back pockets for removable stays. Latex-free elastic. Stays sold separately.

8" back height tapered to a 4" front  
 10" back height tapered to a 6" front  
 12" back height tapered to a 6" front

Optional aluminum stays sold in pairs.

Indications: Mild to moderate low back pain, lumbar sprains and strains, lumbar muscle weakness and lumbar instability.



SIZE	CIRCUMFERENCE
Ex-Small.....	25" – 29"
Small.....	29" – 33"
Medium.....	33" – 37"
Med/Large.....	37" – 41"
Large.....	41" – 45"
Ex-Large.....	45" – 49"
2X-Large.....	49" – 53"
3X-Large.....	53" – 57"